

M-036 Precision Rotation Stage

Piezo Drive Option for Nanometer Precision



M-036.P0 rotation stage with piezo drive

- Sub-Microradian Resolution
- 360° Coarse Range, up to 21° Fine Range
- Precision Micrometer or DC Motor Drives
- Piezo Option for High-Resolution Scanning and Tracking
- Ø 30 mm Clear Aperture

M-036 series precision rotation stages with tangent-arm drive feature high resolution, excellent repeatability and minimum wobble. The stages are equipped with double-row ball bearings for zero backlash and high load capacity. Both the rotation platform and the scale ring (graduated in 2-degree increments) can be independently coarse positioned over 360° degrees and then locked with screws.

Drive Options

A total of six different drive types are offered. They include various combinations of piezoelectric fine-positioners (closed-loop or open-loop), manual and motorized micrometer drives.

Manual Drive

The basic version, the M-036.00, is equipped with a micrometer drive and a zero-

backlash magnetic coupling. The micrometer motion, when converted into rotation, provides a positioning range of 21° degrees (see p. 4-82 for information on how to convert linear input into rotation). The resolution is approximately 15 µrad.

DC Motor Drives

The motorized version, the M-036.D01 features a high-resolution DC motor drive unit (M-227.25, p. 1-42) and has a resolution of 2 µrad. (see p. 7-60 for information on how to convert linear input into rotation). A set of limit switches on the rotation stage protects against overtravel damage.

High-Resolution Piezo Option

For applications requiring extremely high angular resolution, models M-036.PS and M-036.P0 (with manual micrometer drive) and M-036.DS1

and M-036.DP1 (motorized) are available. They have an additional piezoelectric fine adjustment, which can also be used for dynamic operation. The piezo drive has a linear travel range of 45 µm with sub-nanometer resolution, which converts to a rotation range of approx. 0.7 mrad and sub-µrad resolution.

The piezo drives in the M-036.PS and M-036.DS1 versions are also equipped with a position sensor, making closed-loop operation possible, with higher stability, reproducibility and accuracy. For more details on the piezo drives, see the "Piezo Actuators & Components" section (p. 1-61 ff).

Upgrades


M-036 stages without piezo or DC-motor drives can be upgraded at a later date.



Notes

See "Accessories", page 4-89 ff. for adapters, brackets, etc.

Ordering Information

- M-036.00**
Precision Rotation Stage, Ø 100 mm, Micrometer Drive
 - M-036.P0**
Precision Rotation Stage, Ø 100 mm Micrometer + Piezo Drive
 - M-036.PS**
Precision Rotation Stage, Ø 100 mm, Micrometer + Closed-Loop Piezo Drive
 - M-036.D01**
Precision Rotation Stage, Ø 100 mm, DC Motor Drive
 - M-036.DP1**
Precision Rotation Stage, Ø 100 mm, DC Motor + Piezo Drive
 - M-036.DS1**
Precision Rotation Stage, Ø 100 mm, DC Motor + Closed-Loop Piezo Drive
- Upgrade Kits**
- M-036.U0**
Upgrade Kit with Open-loop Piezo Drive
 - M-036.US**
Upgrade Kit with Closed-Loop Piezo Drive
 - M-036.UD**
Upgrade Kit with DC Motor Drive (for factory installation)
- Ask about custom designs!**





Technical Data

Model	M-036.00	M-036.P0	M-036.PS	M-036.D01	M-036.DP1	M-036.DS1	Units
Coarse rotation range	360	360	360	360	360	360	°
Rotation range (micrometer drive)	21	21	21	19**	19**	19**	°
Rotation range (piezo drive)	-	700	700	-	700	700	μrad
Minimum incremental motion (piezo drive)	-	<1	<1	-	<1	<1	μrad
Repeatability (piezo drive)	-	-	2	-	-	2	μrad
Unidirectional repeatability (motor drive)	-	-	-	10	10	10	μrad
Backlash (motor drive)	-	-	-	40	40	40	μrad
Design resolution (motor drive)	-	-	-	0.08	0.08	0.08	μrad
Minimum incremental motion (motor drive)	-	-	-	2	2	2	μrad
Minimum incremental motion (micrometer drive)	23	23	23	-	-	-	μrad
Rotation / linear input	15	15	15	15	15	15	μrad/μm
Tangent-arm length	66	66	66	66	66	66	mm
Wobble	<75	<75	<75	<75	<75	<75	μrad
Max. velocity	-	-	-	0.8	0.8	0.8	°/s
Max. axial force	±400	±400	±400	±400	±400	±400	N
Max. torque (θ _x , θ _y)	±6	±6	±6	±6	±6	±6	Nm
Max. torque CW*	2.6	2.6	2.6	2.6	2.6	2.6	Nm
Max. torque CCW*	0.075	0.075	0.075	0.075	0.075	0.075	Nm
Drive (manual or motor)	M-624	M-624	M-624	M-227.25	M-227.25	M-227.25	
Piezo drive	-	P-840.30	P-841.30	-	P-840.30	P-841.30	
Mass	0.85	0.95	0.97	1.05	1.15	1.17	kg
Body material	Al, St	Al, St	Al, St	Al, St	Al, St	Al, St	
Recommended controllers	-	-	-	C-843, C-848, C-863	C-843, C-848, C-863	C-843, C-848, C-863 (p. 4-120, p. 4-122, p. 4-114)	
Recommended piezo controllers	-	E-660, E-610 E-500 System	E-610 E-500 System	-	E-660, E-610 E-500 System	E-610 (p. 2-110) E-500 System (p. 2-142)	

*CW: clockwise CCW: counter-clockwise

**Limited by limit switch position.

Linear Actuators & Motors

Nanopositioning / Piezoelectrics

Nanometrology

Micropositioning

Hexapod 6-Axis Systems / Parallel Kinematics

Linear Stages

Translation (X)

Vertical (Y)

Multi-Axis

Rotary & Tilt Stages

Accessories

Servo & Stepper Motor Controllers

Single-Channel

Hybrid

Multi-Channel

Micropositioning Fundamentals

Index