

## N-725 PIFOC Piezo Nanofocusing Z Drive with NEXACT® Linear Motor



## Contents

<b>About this Document</b>	<b>3</b>
Symbols and Typographic Conventions.....	3
Other Applicable Documents .....	3
Downloading Manuals.....	4
<b>Safety</b>	<b>4</b>
Intended Use .....	4
Safety Precautions.....	4
<b>Product Description</b>	<b>5</b>
Product View .....	5
Scope of Delivery.....	6
Suitable Electronics .....	6
Accessories .....	6
<b>Unpacking</b>	<b>7</b>
<b>Installation</b>	<b>7</b>
Connecting the N-725 to the Protective Earth Conductor .....	7
Fastening the N-725 .....	8
<b>Startup and Operation</b>	<b>8</b>
General Notes on Startup and Operation .....	8
Starting Up and Operating the N-725.....	9
<b>Maintenance</b>	<b>9</b>
Cleaning the N-725 .....	10
<b>Customer Service</b>	<b>10</b>
<b>Technical Data</b>	<b>11</b>
Specifications.....	11
Maximum Ratings.....	12
Ambient Conditions and Classifications .....	12
Dimensions .....	13
<b>Old Equipment Disposal</b>	<b>14</b>

## About this Document

This user manual contains information necessary for using the N-725 as intended.

It assumes that the reader has a fundamental understanding of basic servo systems as well as motion control concepts and applicable safety procedures.

## Symbols and Typographic Conventions

The following symbols and typographic conventions are used in this user manual:

### NOTICE



#### Dangerous situation

If not avoided, the dangerous situation will result in damage to the equipment.

- Actions to take to avoid the situation.

### INFORMATION

Information for easier handling, tricks, tips, etc.

The following symbols and markings are used in the user manuals of PI:

Symbol	Meaning
1.	Action consisting of several steps with strict sequential order
2.	
➤	Action consisting of one or more steps without relevant sequential order
▪	Bullet
p. 5	Cross-reference to page 5

## Other Applicable Documents

The devices and software tools mentioned in this technical note are described in separate manuals.

Product	Document
P-721.0xQ	P721T0002 technical note
E-861.1A1	PZ205E user manual
PIMikroMove	SM148E software manual

## Downloading Manuals

### INFORMATION

If a manual is missing or problems occur with downloading:

- Contact our customer service department (p. 10).

1. Open the website [www.pi.ws](http://www.pi.ws).
2. Search the website for the product number (e.g., P-725) or the product family (e.g., PIFOC).
3. Click the corresponding product to open the product detail page.
4. Click the **Downloads** tab.

The manuals are shown under **Documentation**. Software manuals are shown under **General Software Documentation**.

5. Click the desired manual and fill out the inquiry form.

The download link will then be sent to the email address entered.

## Safety


### Intended Use

The N-725 is a laboratory device as defined by DIN EN 61010-1. It is intended for indoor use and use in an environment that is free of dirt, oil and lubricants.

According to its design, the N-725 is intended for maximum pull forces of 10 N and maximum push forces of 10 N. The N-725 can be mounted horizontally or vertically.

The N-725 can only be used as intended when completely mounted and connected and only in conjunction with suitable electronics. The N-725 can be operated using the E-861.1A1 NEXACT® servo controller.

### Safety Precautions



## NOTICE



### **Damage due to manual displacement and unallowable forces on the movable part!**

Manual displacement of the movable part (see figure above) and unallowable forces can cause irreparable damage to the piezo modules in the NEXACT® linear motor.

- Do not displace the movable part of the N-725 manually!
- Do not apply a force  $\geq 10$  N on the movable part!

## NOTICE



### **Fixed operating frequency for permanent operation!**

- For permanent operation of the drives do not exceed operating frequencies of 800 Hz.

## NOTICE



### **Unsuitable cables!**

Unsuitable cables can damage the electronics and the N-725.

- Only use cables from PI for connecting the N-725 to the electronics.

## INFORMATION

Extended cables can affect the performance of the N-725.

- If you need longer cables, contact our customer service department.

## Product Description

### Product View




Figure 1: N-725.2A

- 1 Cable exit for piezo voltage
  - 2 Cable exit for sensor signals
  - 3 Housing
  - 4 M4 hole for connecting the protective earth conductor
  - 5 Movable part of the Z drive with threaded aperture
- Z Positive direction of motion

## Scope of Delivery

Product number	Description
N-725.2A	PIFOC piezo nanofocusing Z drive with NEXACT® linear motor, 2 mm travel range, with linear encoder
000036450	M4 screw set for protective earth, consisting of: <ul style="list-style-type: none"><li>▪ M4x8 flat-head screw with cross recess, ISO 7045</li><li>▪ 2 lock washers</li><li>▪ 2 flat washers</li></ul>
P721T0002	Technical note for QuickLock thread option
N725T0002	User manual for N-725.2A (this document)

## Suitable Electronics

Product number	Description
E-861.1A1	NEXACT® controller, 1 channel, with linear encoder

- To order, contact our customer service department (p. 10).

## Accessories

### P-721.xxQ QuickLock thread adapter set

Product no.	Description
P-721.02Q	QuickLock thread adapter set, M26 × 0.75, for objective and microscope
P-721.03Q	QuickLock thread adapter set, M27 × 0.75, for objective and microscope
P-721.04Q	QuickLock thread adapter set, M28 × 0.75, for objective and microscope
P-721.05Q	QuickLock thread adapter set, M32 × 0.75, for objective and microscope
P-721.06Q	QuickLock thread adapter set, M26 × 1/36", for objective and microscope
P-721.08Q	QuickLock thread adapter set, M19 × 0.75, for objective and microscope
P-721.11Q	QuickLock thread adapter set, M25 × 0.75, for objective and microscope
P-721.12Q	QuickLock thread adapter set, W0.8 × 1/36", for objective and microscope

- To order, contact our customer service department (p. 10).

## Unpacking

1. Unpack the N-725 with care.
2. Compare the contents with the items listed in the contract and the packing list.
3. Inspect the contents for signs of damage. If parts are missing or you notice signs of damage, contact PI immediately.
4. Keep all packaging materials in case the product needs to be returned.

## Installation

### Connecting the N-725 to the Protective Earth Conductor


#### INFORMATION

- Observe the applicable standards for connecting the protective earth conductor.

#### INFORMATION

When a Z drive is grounded via its protective earth connection as well as by the shield of the connecting cable for the electronics, ground loops can occur.

- If a ground loop occurs, contact our customer service department (p. 10).

On the N-725, there is an M4 hole for connecting the protective earth conductor. This hole is marked with the symbol for the protective earth conductor . See "Dimensions" (p. 13) for the location of the hole.




Figure 2: Connecting the protective earth conductor (profile view)

- 1 Housing of the N-725
- 2 Flat washer
- 3 Lock washer
- 4 Screw
- 5 Cable lug
- 6 Protective earth conductor

## Requirements

- ✓ You have read and understood the safety precautions (p. 4).
- ✓ The N-725 is **not** connected to the electronics.

## Tools and accessories

- Suitable protective earth conductor: Cross-sectional area of the cable  $\geq 0.75 \text{ mm}^2$
- M4 screw set supplied for the connecting the protective earth conductor (p. 6)
- Suitable screwdriver

## Connecting the N-725 to the protective earth conductor

1. If necessary, attach a suitable cable lug to the protective earth conductor.
2. Use the M4 screw (together with the flat and lock washers) to fix the cable lug of the protective earth conductor to the protective earth connector of the N-725 as shown in the profile view.
3. Tighten the M4 screw with a torque of 1.2 Nm to 1.5 Nm.
4. Make sure that the contact resistance at all connection points relevant for connecting the protective earth conductor is  $< 0.1 \Omega$  at 25 A.

## Fastening the N-725

### Requirements

- ✓ You have read and understood the safety precautions (p. 4).
- ✓ The N-725 is **not** connected to the electronics.

### Fastening the N-725

- Fasten the objective to the thread intended for this purpose.

For detailed instructions refer to the separate technical note P721T0002 for PIFOC QuickLock thread options.

## Startup and Operation

### General Notes on Startup and Operation

#### NOTICE



#### Uncontrolled oscillation!

Oscillations can cause irreparable damage to the Z drive. Oscillations are indicated by a humming and can result from the following causes:

- The load and/or dynamics of operation differ too much from the calibration settings.
- The Z drive is operated near its resonant frequency.
- If you notice oscillations, stop the Z drive immediately.



## INFORMATION

Systems are preconfigured at the factory to achieve optimal positioning accuracy. If you ordered the N-725 as system with controller (e.g., ND72Z2LAQ), suitable parameter sets are already stored in the controller, and the sensor zero point is already adjusted.

- In this case, **do not** assign a stage type in the PC software **and do not** correct the sensor's zero point.  
Otherwise the parameter sets preconfigured at the factory are **lost** and the optimal performance of the system is not provided anymore.

## Starting Up and Operating the N-725

### Requirements

- ✓ You have read and understood the safety precautions (p. 4) and the general notes on startup and operation.
- ✓ You have installed (p. 6) the N-725 correctly and connected it to the electronics.
- ✓ You have read and understood the user manual for the controller (p. 6) used.

### Starting up and operating the N-725

- Follow the instructions in the manual for the controller for startup and operation of the N-725.

## Maintenance

### NOTICE



#### Misalignment from loosening screws!

The N-725 is maintenance-free and precision adjusted.

- Loosen screws only when instructed in this user manual.
- Do **not** open the N-725.

### NOTICE



#### Damage from ultrasonic cleaning!

Ultrasonic cleaning can damage the N-725.

- Do **not** do any ultrasonic cleaning.

---

## Cleaning the N-725

### Requirements

- ✓ You have disconnected the N-725 from the electronics.

### Cleaning the N-725

- Clean the surfaces of the N-725 with a cloth that is dampened with a mild cleanser or disinfectant (e.g., isopropyl alcohol).

## Customer Service

For inquiries and orders, contact your PI sales engineer or send us an email ([service@pi.de](mailto:service@pi.de)).

- If you have any questions concerning your system, provide the following information:
  - Product and serial numbers of all products in the system
  - Firmware version of the controller (if applicable)
  - Version of the driver or the software (if applicable)
  - Operating system on the PC (if applicable)
- If possible: Take photographs or make videos of your system that can be sent to our customer service department if requested.

The latest versions of the user manuals are available for download (p. 4) on our website.

## Technical Data

### Specifications




	N-725.2A	Unit	Tolerance
Active axes	Z		
<b>Motion and positioning</b>			
Integrated sensor	Optical linear encoder		
Closed-loop travel	2000	µm	
Min. incremental motion, closed-loop	5	nm	typ.
<b>Mechanical properties</b>			
Step-and-settle time for a 3 µm step at 200 g payload, 100 nm settling band	<20	ms	
Recommended load*	700	g	max.
<b>Drive properties</b>			
Piezoceramics	NEXACT®		
<b>Miscellaneous</b>			
Operating temperature range	0 to 50	°C	
Material	Aluminum		
Mass	290	g	±5 %
Cable length	1.5	m	±10 mm

\* For dynamic operation. Higher dynamics are possible with a reduced load.

All specifications based on room temperature (22 °C ±3 °C).

## Maximum Ratings

The N-725 is designed for the following operating data:

Maximum operating voltage	Maximum operating frequency	Maximum power consumption
		
45 V	1500 Hz	20 W

## Ambient Conditions and Classifications

The following ambient conditions and classifications for the N-725 must be observed:

Area of application	For indoor use only
Maximum altitude	2000 m
Air pressure	1100 hPa to 0.1 hPa
Relative humidity	Highest relative humidity 80 % for temperatures up to 31 °C Decreasing linearly to 50 % relative humidity at 40 °C
Storage temperature	0 °C to 70 °C
Transport temperature	-20 °C to 70 °C
Overvoltage category	II
Protection class	I
Degree of pollution	1
Degree of protection according to IEC 60529	IP20

## Dimensions

Dimensions in mm. Note that the decimal places are separated by a comma in the drawings.




Figure 3: Dimensions of the N-725 (please order the adapter set separately)

## Old Equipment Disposal

In accordance with EU law, electrical and electronic equipment may not be disposed of in EU member states via the municipal residual waste.

Dispose of your old equipment according to international, national, and local rules and regulations.

In order to fulfil its responsibility as the product manufacturer, Physik Instrumente (PI) GmbH & Co. KG undertakes environmentally correct disposal of all old PI equipment made available on the market after 13 August 2005 without charge.

Any old PI equipment can be sent free of charge to the following address:

Physik Instrumente (PI) GmbH & Co. KG  
Auf der Roemerstr. 1  
D-76228 Karlsruhe, Germany

