


## F-603 Accessories

### Fiber, Objective and Waveguide Holders



Rear (left to right) Microscope Objective Holders: F-603.14, F-603.12 (mounted on P-611.3SF NanoCube® NanoAlignment System) F-603.11, F-603.13

- Mount on a Variety of PI Alignment Systems
- Precision Machined from High-Strength Aluminum/Brass

The F-603 series component holders for fiber optic applications can be combined with the following micro- and nanopositioning systems from PI:

- F-206 6D-Hexapod-Alignment System
- F-130 and F-131 Alignment Systems
- P-611.3OF and P-611.3SF NanoCube® NanoAligners
- A variety of combinations of NanoCubes® with other micropositioning systems

All adapters are equipped with a 3-mm-wide tongue that fits in a slot machined into the platform of the positioning system.

M3 capscrews and miniature cleats supplied with each F-603 adapter are used to fasten these accessories quickly.

#### Ordering Information

- F-603.11**  
Microscope Objective Holder  
W0.8 x 1/36"
  - F-603.12**  
Microscope Objective Holder  
M19 x 0.75
  - F-603.13**  
Microscope Objective Holder  
M25 x 0.75
  - F-603.14**  
Microscope Objective Holder  
M26 x 0.75
  - F-603.20**  
Vacuum Waveguide Mount
  - F-603.21**  
Fiber Holder for FC-Connector
  - F-603.22**  
Ferrule Holder
  - F-603.60**  
V-Groove Fiber Holder with  
Magnetic Clamp
- Ask about custom designs!**

